

Analyses of the oxidative stability of deep-frying oils

Introduction

The quality of edible oils gains more and more in importance. The project presented compares the oxidative stability of conventional **sunflower oil (SF)** with **high-oleic sunflower oil (HOSF)**. HOSF became more interesting for frying applications over the past years.

Material and methods

Oil samples with and without spiked antioxidants (50 mg/kg and 500 mg/kg mixed tocopherols, respectively) were stored at 80°C for 9 days.

The following analyses were performed:

- Measurement of free fatty acids and peroxide value according to §64 LFGB-methods L1300.5 and L1300.6
- Measurement of p-anisidine value with the FOODLABfat device (cdR srl, Florence/Italy)
- Measurement of the antioxidative capacity by TEAC-assay according to Re et al. (1999)
- Measurement of alpha- and delta-tocopherol content by HPLC according to the §64 LFGB-method; eluent was hexane/2-propanol (99/1 v/v)
- Measurement of volatile compounds by headspace-SPME-GC/MS with the following conditions [Jelen et al.(2000)]:

SPME-parameters	
SPME-phase	CAP/DVB/PDMS
enrichment time	90 minutes
enrichment temperature	25°C
desorption time	5 minutes
desorption temperature	260°C

Results

Free fatty acids and peroxide value: Both parameters showed no differences between the oils.

p-Anisidine value: HOSF with added mixed tocopherols showed significant lower values in comparison to SF.

TEAC-assay: Both oils displayed similar results. Addition of mixed tocopherols showed an increasing antioxidative capacity.

Tocopherol content: In the early stage of the accelerated storage test HOSF and SF showed similar contents of alpha-tocopherol. Delta-tocopherol could only be found in supplemented oils. After storage alpha-tocopherol could be detected only in HOSF.

Volatile compounds: 10 major volatile lipid oxidation compounds were selected by their quantitative appearance and their low odour threshold. In both oils hexanal was the only major volatile compound that could be detected in the untreated oils. (E)-2-heptenal was found in considerably higher amounts in SF compared to HOSF. In SF the amount of (E)-2-heptenal decreased with increased tocopherol content. Octanal and nonanal could be found in significant larger amounts in HOSF.




Fig.1: p-Anisidine value of the sunflower oils




Fig.2: Antioxidative capacity of the sunflower oils (trolox equivalents, TEAC-assay)




Fig.3: alpha-tocopherol content in sunflower oils; n.d.: not detectable (<5,61 mg/kg oil)




Fig.4: delta-tocopherol content in sunflower oils; n.d.: not detectable (<7,67 mg/kg oil)




Fig.5: Major volatile compounds in SF oil (above) and HOSF (below)

Tab.1: Volatile compounds in SF oil and HOSF oil (HS-SPME-GC/MS analysis)

Substances	retention index	Substances	retention index	Substances	retention index	Substances	retention index
Alcohol		(Z)-2-Octenal	1049	(E)-2-Hexene	622	Dihydro-2(3H)-furanone	913
Ethanol	450	(E)-2-Octenal	1060	1-Heptene	690	Dihydro-5-methyl-2(3H)-furanone	955
1-Butanol	657	Nonanal	1109	(E)-2-Heptene	706	2-Pentylfuran	992
1-Penten-3-ol	679	(E)-2,4-Decadienal	1117	(E)-3-Heptene	711	Dihydro-5-ethyl-2(3H)-furanone	1055
1-Pentanol	769	(E)-2-Nonenal	1177	1-Octene	799	Dihydro-5-propyl-2(3H)-furanone	1171
1-Hexanol	871	Decanal	1221	(Z)-2-Octene	807		
1-Heptanol	972	(E,E)-2,4-Nonadienal	1230	(E)-2-Octene	815	Ketone	
1-Octen-3-ol	981	(Z)-2-Decenal	1262	Acetic acid	602	Acetone	500
1-Octanol	1073	(E)-2-Decenal	1273	Propanoic acid	691	2-Pentanone	686
Aldehyde		(E)-2,4-Decadienal	1298	Butanoic acid	754	4-Hydroxy-4-methyl-2-pentanone	843
2-Propenal	496	(E)-2,4-Decadienal	1324	Pentanoic acid	895	2-Heptanone	892
Propenal	500	(E)-2-Undecenal	1372	Hexanoic acid	1002	1-Octen-3-one	980
(E)-2-Butenal	646			Heptanoic acid	1081	3-Octanone	988
Pentanal	698	Alkane		Octanoic acid	1189	(E)-3-Octen-2-one	1041
(E)-2-Pentenal	747	Pentane	500	Nonanoic acid	1277	2-Nonanone	1094
Hexanal	800	Hexane	600			(E)-3-Nonen-2-one	1156
(E)-2-Hexenal	854	Heptane	700	Carboxylic Acid		2-Decanone	1207
Heptanal	901	Octane	800	Acetic acid	602		
(E)-2-Heptenal	957	Nonane	900	Propanoic acid	691	Nitrite	
(Z,E)-2,4-Heptadienal	997	Decane	1000	Butanoic acid	754	3-Butenenitrile	654
Octanal	1003	Alkene		Pentanoic acid	895	Sulfurous Substance	
(E,E)-2,4-Heptadienal	1011	1-Hexene	596	Hexanoic acid	1002	Carbon disulfide	537
		(Z)-2-Hexene	607	Heptanoic acid	1081	Terpene	
				Octanoic acid	1189	Limonene	1029
				Nonanoic acid	1277		
				Aromate			
				Toluene	759		
				Ester			
				Ethyl acetate	614		
				Furan			
				2-Ethylfuran	702		